TissuePatch™ Sealant Film

Surgical Sealant and Barrier

Kirurgisk forsegling og barriere

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GB

TissuePatch™ – **Surgical Sealant and Barrier**

Multilaminata adhesive, sealant and barrier film for adjunctive prevention of air, blood and/or fluid leakage

Instructions for Handling and Use

DESCRIPTION

TissuePatch is a STERILE, self-adhesive, absorbable surgical sealant and barrier. The product is a clear film with a blue **TissuePatch** logo visible through the patch. TissuePatch is available in five different sizes summarised below and one size is also available as a multipack.

Size	Ref. No.
50 × 25 × 0,04mm	TP-01
50 × 50 × 0,04 mm	TP-02
50 × 100 × 0,04 mm	TP-03
100 × 100 × 0,04 mm	TP-04
30 × 80 × 0,04 mm	TP-05
30 × 80 × 0,04 mm	TP-06 (4 patch multipack)

For each product size, a single device is packaged within a sterile inner pouch, which is further sheathed in an outer pouch the outside of which is NON-STERILE. TissuePatch is for single patient use only and once opened may not be resterilised. For product should not re-used as this presents numerous hazards including the risk of cross contamination and lack of sterility.

COMPOSITION

TissuePatch is a four layer device comprising alternate layers of poly(lactide-co-glycolide) and poly(N-vinyl-pyrrolidone, -co-acrylic acid, -co-N-hydroxysuccinimide ester of acrylic acid). TissueBond™ TissuePatch contains 0,015% w/w Methylene Blue.

STORAGE

TissuePatch is to be stored at 2-25°C. The product should not be frozen. The use by date of this product is printed on the packaging.

MODE OF ACTION

TissuePatch is a multilayered product. The poly(lactide-co-glycolide) films act as a cohesive matrix within the patch and as an anti-adhesion barrier film on the patch's top surface. The film of TissueBond™ becomes adhesive when applied to tissue and immediately adheres to the surface by ionic and covalent interactions and acts as a sealant to prevent leakage. The patch remains in position while it slowly degrades until substantially absorbed.

INDICATIONS

TissuePatch is indicated for use to seal and reinforce against air leakage in thoracic surgery (open and minimally invasive) and against leakage of low pressure or oozing bleeding or fluid leakage following surgical procedures on soft tissue. TissuePatch is intended for use as an adjunct and is not intended to replace sutures or staples, as appropriate, in tissue approximation.

KONTRA-INDICATIONS

The use of TissuePatch must be considered on an individual basis. It should not be used on patients with a known allergy or sensitisation to its constituent components.

NOT for use to seal defects in neurological applications.

NOT for intravascular use.

WARNINGS AND SAFETY INFORMATION

TissuePatch should not be used as the primary means of support. The product is to be used as an adjunct, providing temporary support.

TissuePatch is not a thrombotogenic haemostat. It should not be used to control arterial bleeds; it can however be used adjunctively to control low pressure/oozing bleeding.

TissuePatch should not be used on mucosal surfaces.

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Send al korrespondance til Deres lokale distributor

DO NOT use TissuePatch if the temperature indicator dot is black; this indicates that the product has been exposed to excessive temperatures and is not suitable for use.

CLINICAL INSTRUCTIONS

The TissuePatch outer package is non-sterile and the wet wrap should be opened by a non-scrubbed person. Care needs to be taken to expose the inner sterile pouch to the sterile surgical field. It is recommended that the inner pouch should be opened and product dispensed immediately prior to use in order to minimize handling and avoid damage (including cutting) to the product. The inner pouch is then opened and the sterile product removed. Care needs to be taken to ensure that the product remains dry before application.

TissuePatch is designed to be an appropriate size for use, but may be cut to size. It is recommended that to avoid damage, TissuePatch should be cut to the appropriate size within the inner pouch. The product should cover a greater area than the area to be treated. For example, during treatment of organs such as liver, spleen or lung allow at least a 1cm margin between the injury and the edge of the patch.

Blood or other fluid leakage should be stemmed to a minimum. Just prior to application of the product to the target tissue, the site should be wiped gently with a dry swab to remove any contaminating materials, e.g. blood and fluid from the tissue.

On application, the blue **TissuePatch** logo should face the user. During application, hold a swab (which has been slightly moistened with saline) evenly over the entire area of the patch for a period of at least 60 seconds to ensure good contact between TissuePatch and tissue surface. TissuePatch is flexible and may be folded and pressed into position.

If following application TissuePatch has not adhered satisfactorily to the target tissue the application described above may be repeated for a further 60 seconds. If the product continues not to adhere it should be removed from the site and the tissue should be prepared as described above for reapplication of a new patch.

If TissuePatch is applied to the incorrect site it may be removed from the site. If further TissuePatch is required, additional product may be applied. Preclinical data supports safe use of up to 200cm² of product.

Size	Maximum Dosage
50 × 25 × 0,04mm	16 = 200cm²
50 × 50 × 0,04mm	8 = 200cm²
30 × 80 × 0,04mm	8 = 200cm²
50 × 100 × 0,04mm	4 = 200cm²
100 × 100 × 0,04mm	2 = 200cm²

Direct all correspondence to your local distributor

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