TissuePatchDural™

sealant film

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Tel. +44 (0)113 200 0500 Fax +44 (0)113 240 7343

GB

TissuePatchDural™ – Dural Sealant and Barrier

DESCRIPTION
TissuePatchDural is a sterile single-use self-adhesive absorbable surgical sealant and barrier. The film of TissueBond™ becomes adhesive when applied to tissue and immediately adheres to the surface by ionic and covalent interactions and acts as a sealant to prevent leakage. The patch remains in position while it slowly degrades until substantially absorbed.

Ref. No.	Size	Total number that may be used (per patient)
TD-01	50 × 25 × 0.04mm	8
TD-02	50 × 50 × 0.04mm	4
TD-03	50 × 100 × 0.04 mm	2
TD-04	100 × 100 × 0.04 mm	1
TD-07	25 × 100 × 0.04 mm	4

MODE OF ACTION
TissuePatchDural is a multilayered product. The poly(lactide-co-glycolide) films act as a cohesive matrix within the patch and as an anti-adhesion barrier film on the patch's top surface. The film of TissueBond™ becomes adhesive when applied to tissue and immediately adheres to the surface by ionic and covalent interactions and acts as a sealant to prevent leakage. The patch remains in position while it slowly degrades until substantially absorbed.

INDICATIONS
TissuePatchDural is indicated for use to seal and reinforce against fluid (including cerebrospinal fluid (CSF)) and/or blood leakage where repair of the dura mater is required. TissuePatchDural is intended for use as an adjunct.

CONTRA-INDICATIONS
Not intended for intravascular use. Not intended to be used as a dura mater substitute. Not intended to replace sutures, staples or clips, as appropriate, in tissue approximation.

WARNINGS AND SAFETY INFORMATION

- TissuePatchDural should not be used on patients with a known allergy or sensitization to its constituent components. The use of must be considered on an individual basis.
- TissuePatchDural should not be used as the primary means of support. The product is to be used as an adjunct, providing temporary support.
- TissuePatchDural is not a thrombotic haemostat. It should not be used to control arterial bleeds. It can however be used adjunctively to control low pressure oozyng/bleeding.
- TissuePatchDural should not be used on the skin.
- If there is any evidence of damage to TissuePatchDural the product shall NOT be used.
- TissuePatchDural should not be used if an active infection is present at the surgical site.
- TissuePatchDural adhesion can be affected by contact with solutions containing iodine. Such substances should NOT be applied to the site before sealing with TissuePatchDural, though can be used following patch application. TissuePatchDural may swell slightly to total thickness of 0.5mm. Care should therefore be taken when using the product in confined bony structures where nerves are present.
- TissuePatchDural may cause a local, mild inflammatory reaction leading to encapsulation by inflammatory cells with some fibrosis.

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Manufactured by Tissueumed Limited

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Hersteller: Tissueumed Limited