TissuePatch™

sealant film

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GB

TissuePatch™

Multilaminatė adhesivė, sealant and barrier film for adjunctive prevention of air, blood and/or fluid leakage

DESCRIPTION
TissuePatch is a sterile single-use self-adhesive absorbable surgical sealant and barrier. The product is a clear, four layer film, comprising alternate layers of poly(lactide-co-glycolide) and poly(N-vinyl-pyrrolidone50-co-acrylic acid25-co-N-hydroxy succinimide ester of acrylate25) ‘TissueBond™’. The blue ‘TissuePatch’ logo contains 50.015% w/v Methylene Blue and is visible through the clear film. A single sterile device is packaged within double barrier pouch packaging.

Ref. No.	Size	Total number that may be used (per patient)	All other applications
TP-01	50 × 25 × 0,04 mm	8	16
TP-02	50 × 50 × 0,04 mm	4	8
TP-03	50 × 100 × 0,04 mm	2	4
TP-04	100 × 100 × 0,04 mm	1	2
TP-05	30 × 80 × 0,04 mm	4	8
TP-06	30 × 80 × 0,04 mm (4 patch multipack)	1	2
TP-07	25 × 100 × 0,04 mm	4	8

MODE OF ACTION
TissuePatch is a multilayered product. The poly(lactide-co-glycolide) films act as a cohesive matrix within the patch and as an anti-adhesion barrier film on the patch's top surface. The film of TissueBond™ becomes adhesive when applied to tissue and immediately adheres to the surface by ionic and covalent interactions and acts as a sealant to prevent leakage. The patch remains in position while it slowly degrades until substantially absorbed.

INDICATIONS
TissuePatch is indicated for use to seal and reinforce against:
• Fluid (including cerebrospinal fluid (CSF)) and/or blood leakage where repair of the dura mater is required;
• Air leakage in thoracic surgery (open and minimally invasive);
• Leakage of low pressure or oozing bleeding or fluid leakage following surgical procedures on soft tissue.
TissuePatch is intended for use as an adjunct and is not intended to replace sutures or staples, as appropriate, in tissue approximation.

CONTRA-INDICATIONS
Not intended for intravascular use.
Not intended to replace sutures, staples or clips, as appropriate, in tissue approximation.
Not intended to be used as a dura mater substitute/replacement.

WARNINGS AND SAFETY INFORMATION
• TissuePatch should not be used on patients with a known allergy or sensitisation to its constituent components. The use of TissuePatch must be considered on an individual basis.
• TissuePatch should not be used as the primary means of support. The product is to be used as an adjunct, providing temporary support.
• TissuePatch is not a thrombogenic haemostat. It should not be used to control arterial bleeds; it can however be used adjunctively to control low pressure oozing/bleeding.
• TissuePatch should not be wrapped around gastrointestinal anastomoses.
• TissuePatch should not be used on the skin.
• If there is any evidence of damage to TissuePatch the product shall NOT be used.
• TissuePatch should not be used if an active infection is present at the surgical site.

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