

TissuePatch™

sealant film

TissuePatch™ is a synthetic, self-adhesive surgical film for protection against leaks




Fast, Safe & Effective Sealing

It Adheres,
It Seals,
It's Absorbed.

Exceptional Strength for Leak Prevention

Performance: Burst Pressure



* Evaluation of Absorbable Surgical Sealants: In vitro Testing. Patrick K. Campbell, PhD; Steven L. Bennett, PhD; Art Driscoll; Amar S. Sawhney, PhD


** Control of bleeding in surgical procedures: critical appraisal of HEMOPATCH (Sealing Hemostat). Med Devices (Auckl). 2015 Dec 22;9:1-10

*** www.tachosil.com/productdetails/efficiency-extensibility-adhesive/carnot-carbondiffuser

¹ Data on file at Tissuemed. Tested in accordance with ASTM F2392.

Clinical Evidence Supporting Reduction in Hospital Stay

TissuePatch was evaluated in a single centre study involving patients undergoing thyroid surgery*



Clinical Publications

Burrows S *et al.* The use of a synthetic, self-adhesive patch for the reinforcement of pharyngeal closure in laryngectomy and pharyngolaryngectomy patients: Initial case series. [Presented at IFOS, Paris, France, June 2017.](#)

*Cheng L *et al.* The use of Tissue Patch, a self-adhesive sealant film to prevent post-operative vascular leakage after thyroid surgery. [Presented at ICOMS, Hong Kong, April 2017.](#)

Surda P *et al.* A purely synthetic and biodegradable material for repair of cerebrospinal fluid rhinorrhoea. [Clin Otolaryngol. 2016 Apr;41\(2\):179-82](#)

Cheng L *et al.* Use of TissuePatch™ sealant film in the management of chyle leak in major neck surgery. [Br J Oral Maxillofac Surg.; 2014 Jan;52\(1\):87-89.](#)

Cheng L *et al.* Thyroid Surgery. [Br J Oral Maxillofac Surg. 2012 Oct;50\(7\):585-591.](#)

Cheng L *et al.* Use of TissuePatch sealant film in reducing vascular leak after thyroid surgery. [Presented at BAOMS, Nice, France, June 2011.](#)

Seals in less than 60 seconds

Ordering Information

TP-01M	TissuePatch 5 x (50 x 25mm)
TP-02M	TissuePatch 5 x (50 x 50mm)
TP-03M	TissuePatch 5 x (50 x 100mm)
TP-04M	TissuePatch 5 x (100 x 100mm)
TP-05M	TissuePatch 5 x (30 x 80mm)
TP-07M	TissuePatch 5 x (25 x 100mm)
Box of 5, Single Packed Patches	